
# INSTRUCTIONS FOR USE: LEIBINGER® MALLEABLE BURR HOLE COVERS

**MANUFACTURED AND DISTRIBUTED BY:**  
**Stryker Leibinger GmbH & Co. KG**  
Bötzing Straße 41  
D-79111 Freiburg, Germany  
Phone: (0761) 4512-0  
Fax: (0761) 4512-120

**DISTRIBUTED IN THE U.S. BY:**  
**Stryker® Leibinger®**  
4100 East Milham Avenue  
Kalamazoo, MI 49001, USA  
Toll free: (800) 962-6558  
Phone: (269) 324-5346

**GB**      **ALL OF THESE INSTRUCTIONS FOR USE MUST BE READ CAREFULLY PRIOR TO CLINICAL USE**      **USA**

**INDICATIONS**

Leibinger® malleable burr hole covers are intended for use in cranial reconstruction surgery, specifically to cover burr holes and secure cranial bone flaps.

**POSSIBLE ADVERSE EFFECTS**

In many instances, adverse results may be clinically related rather than implant related.

1. Osteoporosis, inhibited revascularization, bone resorption can cause loosening of the device.
2. Infection at or near the placement site.
3. Nerve damage may occur as a result of the surgical trauma.
4. Metal sensitivity reactions in patients following implantation have rarely been reported. An evaluation as to their significance is pending further clinical evaluations.

**WARNINGS AND PRECAUTIONS**

Responsibility for proper selection of patients, adequate training, experience in the choice and placement of implants and the decision to leave or remove implants postoperatively, rests with the surgeon.

The surgeon should discuss the expectations of surgery inherent in the use of the product with the patient. Particular attention should be given to a discussion postoperatively and the necessity for periodic medical follow-up.

The correct selection of the product is extremely important. The product should be used in the correct anatomic location, consistent with accepted standards for internal fixation. Failure to use the appropriate product for the application may result in a premature clinical failure. Failure to use the proper component to maintain adequate blood supply and provide rigid fixation may result in loosening, bending, or fracturing of the product and/or bone.

Careful handling and storage of the product is required. Scratching or damage to the component can significantly reduce the strength and fatigue resistance of the product.

Once applied, the product should never be reused. Although it may appear undamaged, previous stresses may have created imperfections that could reduce its service life.

The patient should be advised to report any unusual changes of the operated site to his surgeon. The patient should be closely monitored if a change at the fixation site has been detected. The surgeon should evaluate the possibility of subsequent clinical failure, and discuss with the patient the need for any measures deemed necessary to aid healing.

Because titanium hardness increases and ductility decreases due to cold working during the bending process, it is essential to ensure that the desired shape of the implant is reached in as few bends as possible. Excessive bending can lead to postoperative implant failure.

Deformed screw recesses (countersink in plate) signify not only an increased risk of breakage in these areas, but also impair the accurate fit of the screw head to the plate.

The tab is not intended for implantation.

**CLEANING**

- New products must be carefully cleaned before initial sterilization. Trained personnel must perform cleaning (manual and/or machine cleaning, ultrasound treatment, etc.) along with maintenance and mechanical inspection prior to initial sterilization.
- Exact compliance with the equipment manufacturers' user instructions and recommendations for chemical detergents is required.

**STERILIZATION**

- If not expressly specified as sterile, the product is supplied non-sterile.
- Exact compliance is required with the manufacturers' user instructions for sterilizers.
- It is the responsibility of the user facility to make sure that special cleaning and sterilization methods are used for deactivation of specific pathogens and to validate sterilization cycles to account for differences in sterilization chambers, wrapping methods and load configurations.
- All non-sterile products are sterilizable by steam sterilization (autoclaving). For initial sterilization and resterilization the following parameters can be used:

	Individually wrapped implants or implant container	King Combo & systems with similar configuration	CMF Modular & systems with similar configuration
Cycle Type	Prevac	Prevac	Prevac
Cycle Length	6 min.	15 min.	15 min.
Temperature	270°F (132°C)	270°F (132°C)	270°F (132°C)
Drying Time	5 min.	60 min.	30 min.
Steam Pressure	30 PSI (206 kPa)	30 PSI (206 kPa)	30 PSI (206 kPa)

Container is wrapped using the AAMI (Association for the Advancement of Medical Instrumentation) CSR double wrapping technique. This cycle was developed with only one system in the chamber and using the middle shelf.

**USE OF ORIGINAL PRODUCTS**

Implants and instruments are produced and designed to be used together. The use of products from other manufacturers along with Leibinger® products can involve incalculable risks and/or contamination of the material and misalignments of implant to instrument, thereby endangering the patient, user or third parties.

**HANDLING INFORMATION**

1. Malleable burr hole covers are equipped with a tab that must be removed prior to closure of the surgery site. The tab is not intended for implantation.
2. Malleable burr hole covers can be manually adapted to individual surface contours without the use of bending instruments.

**PLACEMENT AND FIXATION PROCEDURE**

1. Holding the burr hole cover by the attached tab, it is placed over the previously prepared burr hole, fixed on the bone flap side with one screw and lightly tightened.
2. The burr hole cover is rotated onto the bone flap to facilitate the osteotomy of the flap and the flap is removed.
3. Upon completion of the skull procedure, the flap is replaced, the burr hole cover rotated back over the cut, and the final screws are placed.
4. The tab is left in place until the burr hole cover has been fixated at which time the tab is simply removed by bending the tab from the left to the right.

**Caution:** Federal law in the USA restricts the device to sale by or on the order of a physician or hospital.

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